2026.04.17

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Review and approval

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Establish standard operating procedure (SOP) for HLX49-Ab Size Exclusion Chromatographic method (SEC) to determine the content of monomer, high molecular weight species and low molecular weight species of antibody to ensure standardized operation and normal running of the analysis.

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This method is applicable to the determination of SEC purity and impurity for antibody solution, stability samples and other samples of HLX49-Ab manufactured by Henlius.

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HMWS—High Molecular Weight Species

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Analysts should strictly follow this procedure for experimental operation, recording, data processing and result reporting.

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The reviewer should carefully review the accuracy of the record and data calculation.

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The principal of ASD department or authorized designee is responsible for monitoring the implements of the SOP.

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BAD-GE-014 The Operational Procedure of Chromatographic Column for Physichemical and E&L.

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BAD-EQ-003 Agilent High-Performance Liquid Chromatography Operating Procedure.

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BAD-GE-053 Chromatographic Peak Identification and Integration Procedures for ASD Physical-chemical Group.

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Note 1: reagents and materials with the same function can be used with each other.

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Note 2: The guard column can extend the lifespan of the column without affecting the results, it can be used as needed.

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(Note: instruments and equipment with the same function can be used mutually)

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Size-exclusion chromatography (SEC) is a liquid chromatography technique that separates constituents under test by molecular size. SEC separation principle is molecular sieve effect of gel chromatographic columns. A chromatographic column is filled with hydrophilic bonded silica gel whose surface is full of pores with different sizes. When the substances enter the column, different constituents wherein will be separated by molecular sizes. HMWS, Monomer molecules, and LMWS will be eluted sequentially.

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(Note: Reagents preparation described in 6.4.1.- 6.4.2 can be performed in multiples of all solutions, provided that component contents in each preparation remain unchanged.)

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Weigh 80.0 g of sodium hydroxide to 500 mL of ultrapure water, dissolve and mix well.

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The solution may be stored for 3 months at room temperature.

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Mobile phase (100 mmol/L sodium phosphate with 150 mmol/L sodium chloride, pH 6.8±0.2)

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Weigh 6.1 g sodium dihydrogen phosphate dihydrate, 21.8 g of disodium hydrogen phosphate dodecahydrate and 8.8 g of sodium chloride to 950 mL of ultrapure water, and dissolve. Dilute to 1000 mL with ultrapure water. The pH should be measured at 6.8 (If the pH value is not 6.8, adjust the above solution to pH 6.8±0.2 with 4 mol/L sodium hydroxide). Filter the solution through a membrane filter of 0.45 μm and ultrasonicate for 10 min.

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The solution is preserved in a reagent bottle and may be stored for 1 month at 2-8°C. The mobile phase stored for more than 3 days needs to be filtered through a membrane filter of 0.45 μm before using.

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Single preparation for blank, system suitability samples, in-process control samples, research samples and other samples. Duplicated preparation for batch release and stability samples. (The system suitability sample adopts the current reference sample).

(Note: Mixing can be operated by pipetting aspiration or vortexing gently during preparation, pipetting not more than 10 times, vortexing not more than 30 seconds.)

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Preparation of System Suitability Samples

标题: HLX49-Ab 分子排阻色谱检测方法标准操作规程**Title: Standard Operating Procedure of HLX49-Ab Size Exclusion Chromatography Method**

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Dilute appropriate volume of system suitability sample with mobile phase to 1.0 mg/mL, gentle mixing, single preparation.

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Dilute the excipient buffer according to the preparation of the system suitability samples (6.5.1.).

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This method can detect samples in the range of 0.3-2.0 mg/mL, which means sampling amounts of 7.5-50 µg.

6.5.3.1. 上游中控样品和PD步骤样品的前处理**Preparation of In-process control and PD samples for upstream**

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After the sample is purified according to the Standard Operating Procedure (BAD-TP-248), proceed with sample preparation based on sections 6.5.3.2 to 6.5.3.4.

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If the concentration of sample is more than 1.0 mg/mL, dilute sample with mobile phase to 1.0 mg/mL, mix well, and inject 25 µL.

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If the sample concentration was ≥ 0.3 mg/mL and ≤ 1.0 mg/mL, inject 25 µL.

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If the sample concentration was < 0.3 mg/mL and ≥ 0.15 mg/mL, inject 50 µL.

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If the sample concentration was < 0.15 mg/mL, inject 50 µL, and the result is for information only.

标题: HLX49-Ab 分子排阻色谱检测方法标准操作规程**Title: Standard Operating Procedure of HLX49-Ab Size Exclusion Chromatography Method**

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(Note: If necessary, the concentration of the system suitability sample and the test sample can be rounded to the nearest integer for calculation).

6.5.4. 溶液稳定性 Stability of Prepared Sample Solution

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The blank, system suitability sample and test samples prepared above can be stored at 6°C for 48 h.

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System washing and set-up are performed in accordance with Agilent High-Performance Liquid Chromatography Operating Procedure (BAD-EQ-003).

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VWD Detector: wavelength 280 nm, peak width > 0.1 min (2s response time)

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Equilibrate the system if necessary.

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To equilibrate the system, the blank (6.5.2.) or system suitability sample (6.5.1.) could be injected.

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Sequences of in-process control samples, research samples and other samples.

序列顺序 Sequence Order	类型 Type	样品名称 Sample Name	进样次数 Number of Injection
1	系统适应性 System suitability	空白对照 Blank	1
2		系统适应性样品 System suitability samples	1
3	供试品 (1-N) Test sample 1-N	供试品 1 Test sample 1	1
4		供试品 2 Test sample 2	1
.....			
.....		供试品 N Test sample N	1
.....	系统适应性 System suitability	系统适应性样品 System suitability samples	1

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Note: Start the sequence with one injection of blank and system suitability sample, thereafter one injection of the system suitability sample should follow no more than 10 injections of test samples and the last sample injection.

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Sequences of batch release and stability tests.

标题: HLX49-Ab 分子排阻色谱检测方法标准操作规程
Title: Standard Operating Procedure of HLX49-Ab Size Exclusion Chromatography Method

序列顺序 Sequence Order	类型 Type	样品名称 Sample Name	进样次数 Number of Injection
1	系统适应性 System suitability	空白对照 Blank	1
2		系统适应性样品 System suitability samples	2
3	供试品 (1-N) Test sample 1-N	供试品 1-1 Test sample 1-1	1
4		供试品 1-2 Test sample 1-2	1
.....			
.....		供试品 N-1 Test sample N-1	1
.....		供试品 N-2 Test sample N-2	1
.....	系统适应性 System suitability	系统适应性样品 System suitability samples	1

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Note: Start the sequence with one injection of blank and two injections of system suitability samples, thereafter one injection of the system suitability sample should follow no more than 10 injections of test samples and the last sample injection.

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The column is preserved in accordance with The Operational Procedure of Chromatographic Column for Physicochemical and E&L (BAD-GE-014), referring to the preservation and maintenance method of Nano 5 μm column.

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The system is washed in accordance with Agilent High Performance Liquid Chromatography Operating Procedure (BAD-EQ-003).

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Compare chromatograms of the system suitability sample, test sample, and blank control, adjust integration events to ensure that all HMWS、Monomer and LMWS have been integrated except those peaks presented in the blank. Peak of blank excipient is not integrated.

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*Note: If the above parameters are not applicable, reasonable integration can be carried out based on the actual situation, according to the BAD-GE-053 ASD Chromatographic Peak Identification and Integration Procedures, such as adding Minimum Height and Detected Shoulder etc.

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When the substances enter the column, different constituents wherein will be separated by molecular sizes. HMWS, monomer molecules and LMWS will be eluted sequentially. The largest chromatographic peak is the Monomer component, which is preceded by HMWS and followed by LMWS.

6.11.4. 方法保存 Saving method

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Name and save the processing method.

**标题: HLX49-Ab 分子排阻色谱检测方法标准操作规程****Title: Standard Operating Procedure of HLX49-Ab Size Exclusion Chromatography Method**

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Peak areas% are calculated using area normalization method, while the calculation formulas are detailed as 6.15.

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The sample set method, chromatograms report, processing method and instrument method should be saved and printed for batch release and stability samples, in-process control samples, special research samples.

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The sample set method, chromatograms report should be saved and printed for other samples.

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The blank appearing at the location of target peak must have no interference peak or the interference peak area of blank appearing at the location of target peak is not more than 1.0% of the total peak area of the first injection system suitability sample.

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The theoretical plate count of the monomer peak of system suitability samples in the same sequence must be no less than 4000.

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The RSD of monomer% of all system suitability samples in the same sequence must be no more than 2.0%, and the RD of monomer% must be no more than 2.0% if there are only two-needle system adaptable samples.

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The LOQ of the reduced method is 0.02% and the LOD is 0.01%. If the area percentage

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Title: Standard Operating Procedure of HLX49-Ab Size Exclusion Chromatography Method

is less than 0.02% but not less than 0.01%, reporting '<0.02%' for reference. If the area percentage is less than 0.01%, reporting 'ND'.

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For in-process control samples, research samples and other samples, reporting the HMWS%, Monomer% and LMWS% of each injection. Data is accurate to 0.1%. If 0.0% or 0.00% occur after the result has been modified, the report result is allowed to retain one significant digit.

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For batch release and stability samples, reporting the mean of HMWS%, Monomer% and LMWS%. Data is accurate to 0.1%. If 0.0 or 0.00 occur after the result has been modified, the report result is allowed to retain one significant digit.

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
Figure 1 Typical chromatograms of HLX49-Ab representative sample (Full-scale)

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Figure 2 Typical chromatograms of HLX49-Ab representative sample (Expanded)

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Note: English for reference only, do not fill it.

检验日期: Test date		检验记录编号: Test record number
序号 No.	ASD样品号 ASD Sample Number	样品名称 Sample Name
1		
2		
3		
4		
5		

1. 溶液信息 Solution Information

样品类型 Sample type	厂家/货号 Manufacturer and Cat. No.	批号/配制编号 Batch No. /Preparation No.	有效期至 Expiration Date
系统适应性样品 System Suitability Sample	Henlius		
辅料溶液 Excipient solution	实验室配制或分装 Laboratory preparation or sub package		
流动相 Mobile Phase	实验室配制 Laboratory preparation		

2. 样品制备 Sample Preparation

2.1. 系统适应性样品、空白对照溶液制备 System Suitability Samples and Blank Preparation

样品类型 Sample Type	浓度 Concentration (mg/mL)	样品加入量 Sample Volume (μL)	流动相加入量 Mobile Phase Volume (μL)	制备方式 Preparation Method
空白对照溶液 Blank	N/A			单份制备 Prepared in single
系统适应性样品 System Suitability Sample				单份制备 Prepared in single

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CONFIDENTIAL

2.2. 样品制备 Test Sample Preparation

序号 No.	浓度 Concentration (mg/mL)	样品加入量 Sample Volume (μL)	流动相加入量 Mobile Phase Volume (μL)	制备方式 Preparation Method
1				☐ 单份制备 Prepare in single ☐ 双份制备 Prepare in duplicate
2				☐ 单份制备 Prepare in single ☐ 双份制备 Prepare in duplicate
3				☐ 单份制备 Prepare in single ☐ 双份制备 Prepare in duplicate
4				☐ 单份制备 Prepare in single ☐ 双份制备 Prepare in duplicate
5				☐ 单份制备 Prepare in single ☐ 双份制备 Prepare in duplicate

3. 仪器和色谱柱 Instrument and Column

仪器/色谱柱 Instrument/Column	厂家/型号/参数 Manufacturer/Cat. No./ Parameter			设备编号/序列号 Equipment Number/Serial No.
HPLC	Agilent 1260			
	柱温 Column temperature (°C)	流速 Flow rate (mL/min)	检测波长 Wavelength (nm)	
色谱柱 Colum	Nano BioCore SEC-300 (150×7.8 mm, 3 μm)			

4. 结果 Results

4.1. 结果保存 Save Path

Project Name	
Sample Set Name	
Acq. Method Name	
Processing Method	

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Title: Test Record of HLX49-Ab Size Exclusion Chromatography Method

考察内容 System Suitability Test	结果 Result	结论 Conclusion
空白对照色谱图在目标峰出峰位置无干扰, 或空白对照色谱图于目标峰出峰位置的色谱峰面积不大于第一针系统适应性样品总峰面积的 1.0%。 The blank appearing at the location must have no interference peak or the interference peak area is not greater than 1.0% of the total peak area of the first injection of system suitability sample.	☐ 无干扰峰 no interference peak ☐ 有干扰峰 interference peak, Area%: _____	☐ 通过 Pass ☐ 不通过 Fail

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System Suitability results (Peak area% data accurate to 0.1%. If 0.0% or 0.00% occur after the result has been modified, the report result is allowed to retain one significant digit)

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考察内容 Contents of investigation	结论 Conclusion
同一序列中系统适应性样品单体峰峰面积%的 RSD $\leq$ 2.0%, 如只有两针系统适应性样品, 单体峰面积%的 RD $\leq$ 2.0% The RSD of monomer% of all system suitability samples in the same sequence must be no more than 2.0%, and the RD of monomer% must be no more than 2.0% if there are only two-needle system adaptable samples.	☐ 通过 Pass ☐ 不通过 Fail
每针系统适应性样品单体峰理论塔板数 (EP) $\geq$ 4000 The theoretical plate count (EP) of the monomer peak of system suitability samples in the same sequence must be no less than 4000.	☐ 通过 Pass ☐ 不通过 Fail

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Results: Report data is accurate to 0.1% (Peak area% data accurate to 0.1 %. If 0.0% or 0.00% occur after the result has been modified, the report result is allowed to retain one significant digit).

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Note: The LOQ of the reduced method is 0.02% and the LOD is 0.01%. If the HMWS% and the LMWS% are less than 0.02% but not less than 0.01%, reporting ‘<0.02%’ for reference. If the HMWS% and the LMWS% are less than 0.01%, reporting ‘ND’.

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Sequence table refer to attachment_____, total pages_____.

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















Original chromatograms refer to attachment_____, total pages_____.

6. 备注 Note

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
复核人 Reviewed by: _____ 日期 date: _____

附件编号 Attachment No. PD-GE-014B	版本 Ver. 03	Page 1 of 1
标题：研发部门培训记录 Title: Training Record of R&D Department		

培训主题 Training Subject	HLX49-Ab 分子排阻色谱法检测方法标准操作规程 HLX49-TP-002, Ver.02		
培训人签名/日期 Trainer /Date	 2026.04.17	培训地点 Training Location	松江一厂
培训日期 Training Date	2026.04.17	培训课时 Training Hours	0.5h
培训类型 Training Type	☒ 文件培训 Document training ☐ GxP ☐ 外部培训 External training ☐ 其他 Other: _____		
培训方式 Training Method	☒ 阅读理解 Reading Comprehension ☐ 课堂培训 Classroom Training ☐ 操作培训 Operation Training ☐ 其他 Other: _____		
评估方式* Assessment Type	☐ 书面测试 Written Test ☒ 口头问答 Oral Q&A ☐ 电子测试 Electronic Test ☐ 操作 Operation ☐ 无需评估 No Assessment ☐ 其他 Other: _____		
考核标准 Assessment Criteria	测试正确率 100%通过为合格。		
部门 Department	受训人签名/日期 Trainee's Signature/Date	部门 Department	受训人签名/日期 Trainee's Signature/Date
ASD	 2026.04.17	ASD	 2026.04.17
ASD	 2026.04.17	ASD	 2026.04.17
ASD	 2026.04.17	ASD	 2026.04.17
ASD	 2026.04.17	ASD	 2026.04.17
ASD	 2026.04.17	ASD	 2026.04.17
ASD	 2026.04.17	ASD	 2026.04.17
ASD	 2026.04.17	ASD	 2026.04.17
ASD	 2026.04.17	ASD	 2026.04.17
ASD	 2026.04.17	N/A	N/A

附件编号 Attachment No. PD-GE-014B	版本 Ver. 03	Page 1 of 1
标题：研发部门培训记录 Title: Training Record of R&D Department		

培训是否合格 Pass or not	☒ 是 Yes (以上所有受训人均通过培训 All the above trainees pass the training) ☐ 否 No (请于下方空白处列出未达到考核标准者, 除下方所列人员, 其他受训人均通过培训 Please list in the space below those who do not meet the assessment criteria, except for those listed below, all other trainees pass the training) 评估人/日期:  Evaluated By/Date: 2016.04.17
备注 Notes	N/A

*口头问答将相应的问题和答案附后, 电子测试将学习记录附后, 操作将实验记录或其复印件附后。
 Oral Q&A attaches the questions and answers. Electronic test attaches learning records. Operation attaches experimental records or their copies.